

Informations générales

Piles Zinc-Air A312
Voltage 1,4V

Réglages techniques

		Gain 11 DB	Gain 20 DB
Niveau de sortie	Max	99 dB	99 dB
	HF-Moyenne	93 dB	93 dB
Amplification	Max	15 dB	25 dB
	HF-Moyenne	11 dB	20 dB
Distorsion THD	500Hz	0.5%	0.8%
	1000Hz	0.1%	0.1%
	1600Hz	0.4%	0.3%
Equivalent input noise		37 dB	30 dB
Consommation courant (mA@1,3V)		0.83	0.84

Courbe de gains perçus

